




## *3D Slicer Training Compendium*

---

# Using Plastimatch for Registration and Warping

Tutorial Version: Dec 29, 2015



Gregory C. Sharp  
[gcsharp@partners.org](mailto:gcsharp@partners.org)

Department of Radiation Oncology  
Massachusetts General Hospital





# *Learning Objective*

---

This tutorial is a step-by-step guide, and includes:

- 1) Downloading the Plastimatch extension to 3D Slicer
- 2) Loading the DICOM and DICOM-RT data
- 3) Running and viewing deformable registration
- 4) Exporting DICOM-RT data *(to be developed)*

The 3D Slicer web site is: <http://slicer.org>

The plastimatch web site is: <http://plastimatch.org>

---



# *Prerequisites*

---

Tutorial sample data is required.  
You can get the data from here:


<http://bit.ly/1TnQRv7>

The tutorial “Slicer4 Data Loading and 3D Visualization”  
is recommended, but not required.

<http://bit.ly/1mgOUXh>


## Part 1: Downloading the Plastimatch Extension

Plastimatch is distributed with the SlicerRT extension



Start up 3D Slicer

Choose “Extension Manager”  
from the “View” menu



3D Slicer 4.5.0-1

File Edit View Help



Modules: Welcome to Slicer



S: 0.000mm 1

Extensions Manager

Manage Extensions (0)

Install Extensions

Search...



INSTALL

INSTALL

INSTALL

INSTALL

INSTALL



PETTumorSegmenta..  
Christian Bauer (Univ..  
★★★★★ (0)

SlicerRT  
Csaba Pinter (PerkLa..  
★★★★★ (0)

PercutaneousApproa..  
Koichiro Murakami (S..  
★★★★★ (0)

PetSpectAnalysis  
Martin Bertran, Nata..  
★★★★★ (0)

PickAndPaintExtensi..  
Lucie Macron (Univer..  
★★★★★ (0)


INSTALL

INSTALL

INSTALL

INSTALL

INSTALL



Q3DC

Resection Planner  
Lougheed (Que...  
★★★★★ (0)

PerkTutor  
Tamas Ungi, Matthew..  
★★★★★ (0)

Reporting  
Andrey Fedorov (SPL...  
★★★★★ (0)

INSTALL

INSTALL

INSTALL


INSTALL

INSTALL

Find the SlicerRT extension  
Then click "INSTALL"

Restart

Close



## Extensions Manager

Manage Extensions (1)

Install Extensions

Search...

INSTALL

INSTALL

INSTALL

INSTALL

INSTALL



PETTumorSegmenta..  
Christian Bauer (Univ..  
★★★★★ (0)

INSTALL




SlicerRT  
Csaba Pinter (PerkLa..  
★★★★★ (0)

UNINSTALL




PercutaneousApproa..  
Koichiro Murakami (S..  
★★★★★ (0)

INSTALL




PetSpectAnalysis  
Martin Bertran, Nata..  
★★★★★ (0)

INSTALL



PickAndPaintExtensi..  
Lucie Macron (Univer..  
★★★★★ (0)

INSTALL



Q3DC  
Lucie Macron (Univer..  
★★★★★ (0)

INSTALL



PkModeling  
Yingxuan Zhu (GE), J..  
★★★★★ (0)

INSTALL



ResectionPlanner  
Matt Loughheed (Que..  
★★★★★ (0)

INSTALL



PerkTutor  
Tamas Ungi, Matthew..  
★★★★★ (0)

INSTALL



Reporting  
Andrey Fedorov (SPL..  
★★★★★ (0)

INSTALL


Click "Restart" to restart Slicer

\* Restart requested

Restart

Close








## Part 2: Loading the DICOM and DICOM-RT data



# Welcome

 Load DICOM Data Load Data Customize Slicer Download Sample Data

## Feedback



Share your stories with us and let us know about how 3D Slicer has enabled your research.

We are always interested in improving 3D Slicer, and every submission will be carefully read.

See more at <http://goo.gl/6BvcHm>.

► About

► The Main Window

► Loading and Saving

## Data Probe

Green RAS: ( 291.0, -1.0, -15.4) Coronal Sp: 1.0

L None


F None

B None


Start up 3D Slicer

Click the DCM button






The DICOM Browser window  
will appear. Click "Import".




A file picker dialog will appear.  
Choose the tutorial directory.

Then click "Import"



Another dialog appears.  
Choose "Add Link"



Help & Acknowledgement

Servers

Start Listener

Start Listener when Slicer Starts

DICOM Database and Networking

Show DICOM Browse

Recent DICOM Activity

5 series added to database in the past hour

Today: for MODEL^POPI  
Today: for MODEL^POPI  
Today: for MODEL^POPI  
Today: for MODEL^POPI  
Today: for MODEL^POPI

Refresh

Data Probe

Show Zoomed Slice

L  
F  
B

DICOM Directory Import

Directory import completed.  
2 New Patients  
2 New Studies  
116 New Series  
285 New Instances

OK

Import DICOM files from directory ...

HShome/gcs6/work/plastimatch-tutorial

popi-plan-2015-12-29  
3D\_Slicer\_Plastimatch\_DICOM\_RT\_Tutorial.ppt

Directory: popi-plan-2015-12-29

Files of type: Directories

Copy on import

Import


Cancel

Yet another dialog appears. Click "OK".


10 cm

10 cm









## Part 3: Running and Viewing Deformable Registration

The screenshot shows the 3D Slicer 4.5.0-1 interface. The left sidebar contains the 'Modules' panel, which is currently set to 'DICOM'. A red arrow points to the 'DICOM' module in the list. The main window displays a 3D view of a brain scan with a red arrow pointing to the 'Plastimatch' module in the 'DICOM' list. A second red arrow points to the 'B-spline deformable registration' option in the 'Plastimatch' sub-menu. The bottom right corner shows a 3D view of a brain scan with a red arrow pointing to the 'B-spline deformable registration' option in the 'Plastimatch' sub-menu.

Modules: **DICOM**

- All Modules
- Annotations
- Data
- DataStore
- DICOM
- Editor
- Markups
- Models
- Scene Views
- Segmentations
- Subject Hierarchy
- Transforms
- View Controllers
- Volume Rendering
- Volumes
- Welcome to Slicer
- Radiotherapy
- Plastimatch**
  - DICOM Registration Export
  - Plastimatch
    - Registration**
      - B-spline deformable registration
      - LANDWARP Landmark deformable registration
  - Wizards
  - Informatics
  - Registration
  - Segmentation
  - Quantification
  - Diffusion
  - IGT
  - Filtering
  - Surface Models
  - Converters
  - Endoscopy
  - Utilities
  - Developer Tools
  - Legacy
  - filter

Click the modules menu

Navigate to Plastimatch

"B-spline deformable registration"



3D Slicer 4.5.0-1

File Edit View Help

DATA DCM SAVE Modules: B-spline deformable registration

3DSlicer

Help & Acknowledgement

B-spline deformable registration

Parameter set: B-spline deformable registration

Input/Output

Fixed volume303: Unknown

Moving volume1: Unknown

Input B-spline transform (MRML)None

Input vector field (MRML)None

Input transform (file)None

Fixed fiducialsSelect a MarkupsFiducial

Moving fiducialsSelect a MarkupsFiducial

Output volumeNone

Output B-spline transform (MRML)None

Output vector field (MRML)

Restore Defaults AutoRun

Data Probe

Show Zoomed Slice

L

F

B

R

S: -24.500mm

1

B: 303: Unknown

10 cm

B: 303: Unknown

10 cm

For the tutorial, we will map the dose from image “1:Unknown” onto image “303:Unknown.”

Set fixed to “303”, and moving to “1.”

For “Output B-spline transform (MRML)” option, choose “Create new BsplineTransform”

3D Slicer 4.5.0-1

File Edit View Help

Modules: B-spline deformable registration

3DSlicer

Help & Acknowledgement

B-spline deformable registration

Parameter set: B-spline deformable registration

Input/Output

Fixed volume	303: Unknown
Moving volume	1: Unknown
Input B-spline transform (MRML)	None
Input vector field (MRML)	None
Input transform (file)	None
Fixed fiducials	Select a MarkupsFiducial
Moving fiducials	Select a MarkupsFiducial
Output volume	None
Output B-spline transform (MRML)	Output B-spline transform (MRML)
Output vector field (MRML)	None

Restore Defaults AutoRun

Status: Idle

Click "Apply"

Click "Apply"

3D Slicer 4.5.0-1

File Edit View Help

Modules: B-spline deformable registration

3DSlicer

Help & Acknowledgement

B-spline deformable registration

Parameter set: B-spline deformable registration

Input/Output

Fixed volume: 303: Unknown

Moving volume: 1: Unknown

Input B-spline transform (MRML): None

Input vector field (MRML): None

Input transform (file): None

Fixed fiducials: Select a MarkupsFiducial

Moving fiducials: Select a MarkupsFiducial

Output volume: None

Output B-spline transform (MRML): Output B-spline transform (MRML)

Status: Completed 100%

Restore Defaults AutoRun Cancel Apply

Data Probe

Show Zoomed Slice


L F B

After registration is completed, we want to inspect the results.

The first step is to make the Window/Level settings the same for both images.

Choose the "Volumes" icon





Set the “Active Volume” to “1”

Then click on the “Lung” display setting icon

Repeat this for volume “303” and set it to “Lung” display








Image "1:Unknown" has been transformed.

You may notice the structures don't line up. To transform the structures, click "1:RTSTRUCT: POPI", then click the right arrow.

Next, set the display to  
“Green Only”

Information

Edit

Identity Invert

Display

Apply transform

Transformable:

RTDOSE [1]: Lung  
Dose\_ColorTable  
1: RTSTRUCT: POPI\_ColorTable  
AP a01  
BeamModel\_AP a01  
RL a01  
BeamModel\_RL a01  
LPI  
BeamModel\_LPI  
LPS  
BeamModel\_LPS  
303: Unknown


Transformed:

1: RTSTRUCT: POPI  
1: Unknown

Data Probe

Show Zoomed Slice

L  
F  
B





3D Slicer 4.5.0-1

File Edit View Help

Modules: Transforms

Click the “pin icon”.  
Then click the “>>” button.  
Then click the “Link” button.  
Then set the foreground image to “1:Unknown”

Coronal

1: Unknown  
303: Unknown

1: Unknown

1: Unknown (100%)  
B: 303: Unknown

10 cm

3D Slicer 4.5.0-1

File Edit View Help

Modules: Transforms

Coronal

1.00 1.00 1.00

None

1: Unknown


303: Unknown

A: 0.498mm

You are now ready to begin visualization.

To blend between the two images, drag the blending puck up and down.


Or Ctrl+Left click and drag up-down.



F: 1: Unknown (100%)  
B: 303: Unknown

10 cm

The registration result is fair.  
The diaphragm is aligned, but tumor and  
vessels need improvement.



3D Slicer 4.5.0-1

File Edit View Help

Modules: B-spline deformable registration

3DSlicer

Help & Acknowledgement

Landmark Penalty: 0.005

Max iterations: 50

Output volume\_1: None

Stage 2 Options

Enable Stage 2? ☒

Image subsampling rate (vox): 2,2,1

Grid size (mm): 50.00

Regularization: 0.005

Landmark Penalty: 0.005

Max iterations: 50

Output volume\_2: None

Stage 3 Options

Status: Completed

100%

Restore Defaults AutoRun Cancel Apply

Data Probe

Show Zoomed Slice

L

F

B

Coronal

A: 0.498mm

We will improve the registration by letting it run longer. Return to the B-spline deformable registration module.


Click "Enable Stage 2?"

Then click "Apply"

F: 1: Unknown (0%)  
B: 303: Unknown

10 cm





## Part 4: Exporting DICOM-RT Data

*(This part of the tutorial is under construction)*



# Conclusion

---

Congratulations! You have completed the tutorial.

Please send corrections or suggestions to:

Greg Sharp

[gcsharp@partners.org](mailto:gcsharp@partners.org)

Or visit the web page at:

<http://plastimatch.org>





# *Acknowledgements*

---



**National Alliance for Medical Image Computing**  
NIH U54EB005149

**National Institutes of Health**

NIH / NCI 6-PO1 CA 21239

Federal share of program income earned by MGH on C06CA059267

**Progetto Rocca Foundation**

A collaboration between MIT and Politecnico di Milano

